


# “MAXXCLAVE” - CLASS B AUTOCLAVE



## Unique Features

- Simple and user friendly machine with 4 keys for intuitive navigation.
- Machine as per European B Class Standard with 3 Pulse Pre- Vacuum for effective air removal. Vacuum goes upto -0.90 Bar
- Backlight LCD display to control process.
- Automatic Door Locking for extra safety.
- 8 Pre-programmed Sterilizations cycles & 2 Test cycle (VACUMM/HELIX/Bowie Dick).
- Highly efficient Steam generator to produce high quality saturated steam.
- Automatic water feeding & draining system.
- Excellent post vacuum drying to give dry load in every cycle.
- Inbuilt Thermal printer to maintain records.
- The machine complies with EN 13060, ISO 13485, ISO 9001 & CE which guarantees Safe & Accurate sterilization.



**STERIMAXX INC**

*"A reliable partner in Sterilization !!!"*

# “MAXXCLAVE” - CLASS B AUTOCLAVE

Model	MAXXCLAVE - 17B	MAXXCLAVE - 23B
Chamber Dimensions(mm)(Φ × D)	25X35	25 × 49
Chamber Capacity(liter)	17L	23L
Overall Size(cm)(L × W × H)	48 × 53 × 42	48 × 69 × 42
Heater	1300W	1800W
Net Weight(KGS)	40	44
Gross Weight(KGS)	48	53
Packing Size(cm)	63 × 60 × 50	83 × 60 × 50
Minimum Desk Size(cm)	48 × 40	48 × 50
Power Supply	230V,50/60HZ	
Sterilize Temperature	121°C/134°C	
Sterilize Time	Please refer to different programs	
Dry Procedure	vacuum drying	
Sterilizer Construction	Stainless steel SUS#304	
Water Tank	internal	
Stainless steel instrument containe	3trays	

Independent test program lets you  
you know the machine  
in the optimal condition at any time

Program	Vacuum test	Helix/bdtest
Temperature ( °C )	—	134
Pressure(bar)	-0.8	2.1
Sterilizing time	—	18
Pre-vacuum times	—	3
Vacuum dry	—	10
Total time (min)	8	58



**8** programs and **2** test program  
(Vacuum test, Helix BD test)

Various programs to satisfy your different requirements

Program	Temperature	Pressure	Sterilizing Time	Vacuum Dry	Total Time (min)	Sterilizing Products
porous	134	2.1	4	yes	35-48	Wrapped hollow instruments, porous materials
hollow	134	2.1	4	yes	30-42	Unwrapped hollow instruments
wrapped	134	2.1	4	yes	25-35	Wrapped hollow instruments, porous materials
soild	134	2.1	4	yes	25-35	Unwrapped solid instruments
porous	121	1.1	20	yes	45-55	Wrapped hollow instruments, porous materials
hollow	121	1.1	20	yes	40-50	Unwrapped hollow instruments
wrapped	121	1.1	20	yes	35-48	Unwrapped solid instruments
prion	134	2.1	18	yes	50-60	Unwrapped hollow instruments



**STERIMAXX INC**

"A reliable partner in Sterilization !!!"



A-55, Shankesh Building, Sanghvi Nagar, Near TVS Showroom,  
Mira Bhayander Road, Mira Road (East), Thane 401 107, Maharashtra, India.



+91 998751 4946, 730412 4690



sterimaxx@gmail.com  www.sterimaxx.com