

3M™ Rapid Multi-Enzyme Cleaner

- ✓ Suitable for Manual & ultrasonic cleaning
- ✓ Compatible with Surgical & Medical Instruments
- ✓ Neutral pH
- ✓ Fast Acting Formulation
- ✓ Unique Stable formulation.
- ✓ Contains 4 enzymes - Protease, Amylase, Lipase & Cellulase
Making it highly effective against proteins, carbohydrates, lipids & mucopolysaccharides.



Product Information:

3M™ multi-enzyme cleaners are cleaning solutions, used for cleaning medical and surgical equipment. They have a long shelf life and are environment-friendly. Save time and effort cleaning medical equipment with the 3M™ Rapid Multi-Enzyme Cleaner. It is a fast and easy cleaning solution for surgical and dental instruments, including flexible and rigid endoscopes. This surgical instrument cleaner has a unique stable formulation that keeps the solution active for up to two years from the manufacturing date, when stored at the proposed temperatures between 12 - 25 degrees Celsius. This rapid multi-enzyme cleaner can efficiently digest human secretion and is highly beneficial against proteins, carbohydrates, lipids and mucopolysaccharides. The neutral pH of this multi enzyme cleaner allows it to be easily disposed through a sink. It is suitable for use during ultrasonic bath or manual cleaning. It can be used on all metallic and plastic, medical and dental instruments including flexible and rigid endoscopes, ultrasonic probes, hospital linen and reusable catheters.

STERIMAXX INC

"A reliable partner in Sterilization !!!"

A/55 Shankesh Bldg, Sanghvi Nagar,
Mira Bhayander Road, Mira Road (East).
Thane 401 107 Maharashtra, INDIA.
Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com