

# Water Distiller

- ✓ Easy To use & maintain.
- ✓ Most Economical & 100 % Distill water
- ✓ Runs without running water.
- ✓ In 4.5 Hrs we get 4 Ltrs of Distill Water.
- ✓ Table top mounting.
- ✓ Automatically switches off at the end of cycle.
- ✓ Used for Table top Autoclaves & Ultrasonic Cleaner



## Why Distill Water?

In Autoclave we must use distilled water as the working medium. Even pure water it also contains various minerals, evaporation at high temperatures will produce scale. Over time, the accumulation of scale will be causing the steam leak, causing a small pipe jam, pressure and temperature sensor failures etc.

While in chamber with the whole sterilizing equipment, scale easy to jam into the small pipes and bearings, to make equipment performance degradation, which affects using life. This shows the importance of using distilled water. Clinics, hospitals, laboratories often use a little of distilled water from market which doesn't give any purity and quality assurance even the costing is high at sometime.

Hence, it's the best way to make distilled water by ourselves .You just need to fill into the ordinary tap water, and then you can get high-quality distilled water.

## Technical Specification:

Voltage	AC 220V/110v ±10%
Frequency	50/60Hz±1Hz
Wire plug	750W
Max Hp	4Litres
Tank Capacity	One Ltr/Hr
Class	Φ230×380mm
Sterilizing temperature	Φ180×200mm
Inside Dimension	290×260×450mm
Net weight	5.0kg



**STERIMAXX INC**

*"A reliable partner in Sterilization !!!*

A/55 Shankesh Bldg, Sanghvi Nagar,  
Mira Bhayander Road, Mira Road (East).  
Thane 401 107 Maharashtra, INDIA.  
Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com