

VERTICAL AUTOCLAVE

DOUBLE DRUM DOUBLE JACKETED Automatic Vertical Autoclave

Unique Features:

- ❖ Unique Double wall Design.
- ❖ Heavy Duty Jacket.
- ❖ Fully Automatic Operation.
- ❖ Microprocessor based PID Controller.
- ❖ Dual Display for Temperature & Time.
- ❖ Purge valve Available.
- ❖ Reduced loading height.
- ❖ Inbuilt water level sensor & Temp sensor.



Technical Specification:

Model	BTE 21	BTE 35 FA	BTE75A
Capacity	21 Ltrs	35 Ltrs	75 Ltrs
Chamber Dimension	300 x 300mm (12 x12")	300x500mm (12 x 20 inches)	400x600mm (16 x 24 inches)
Construction	SS304	SS 304	SS304
Power Consumption	2Kw	2Kw	4Kw
Operating Pressure	15PSI	15 Psi	15 Psi
Operating Temperature	121° C	121° C	121° C
Voltage	230v	230V	230V
Display	NO	Digital	Digital
Shipping Dimensions		800 x 570 x 1110 (mm)	900 x 710 x 1215 (mm)



STERIMAXX INC

"A reliable partner in Sterilization !!!

A/55 Shankesh Bldg, Sanghvi Nagar,
Mira Bhayander Road, Mira Road (East).
Thane 401 107 Maharashtra, INDIA.
Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com