

Horizontal Autoclave

- ✓ Easy to read pressure & Temperature.
- ✓ Smooth operation.
- ✓ Cage glass tube with automatic closing safety valve.
- ✓ Heat resistant silicone Door gasket.
- ✓ Low water protection system, to protect the electric heaters.
- ✓ Safety valve for the jacket/Boiler.
- ✓ Hydraulically test to 2.5 times the designed working pressure.
- ✓ Powerful ejector for drying.
- ✓ Printer available (Optional).
- ✓ Full Vacuum thru Vacuum Pump can be done.
- ✓ Insulated with resin bond fiber glass wool covered with M.S / S304 sheet.
- ✓ Self locking safety door cannot be opened while chamber is under pressure.
- ✓ Multi Port operating Valve to control entire cycle operation permits two exhaust speeds, fast & slow.



Model	HCS 80 High Speed	HCS 186 High Speed	HCS 80	HCS 145	HCS 186	HCS 227	HCS 270	HCS 510
Chamber Vol	80	186	80	145	186	227	270	510
Chamber (DIA)	400 (16)	500 (20)	400 (16)	400 (20)	500(20)	500(20)	600(24)	750(30)
mm/inches (DEP)	600 (24)	900 (36)	600 (24)	900 (36)	900(36)	1100(44)	900(36)	1100(44)
Power	18Kw	18Kw	6Kw	9Kw	9Kw	9Kw	9Kw	18Kw
Boiler Cap (in Lit)	45	69	26	44	44	44	44	80
Electrical Req.	42A	40A	24A	24A	24A	24A	24A	63A
Electrical Conn.	3	3	3	3	3	3	3	3
Floor Space WID	1130	1175	1130	1130	1175	1175	1230	1460
Floor Space Dep	1900	2210	1900	2275	2210	2350	2400	3090



STERIMAXX INC

"A reliable partner in Sterilization !!!"

A/55 Shankesh Bldg, Sanghvi Nagar,
Mira Bhayander Road, Mira Road (East).
Thane 401 107 Maharashtra, INDIA.

Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com