

Quick Autoclave

- ✓ Simple Door Locking System (with door interlock)
- ✓ All contact parts are made of SS 304
- ✓ Outer body is made up of high grade plastic
- ✓ Water should be filled from top. Using less than 150 ml.
- ✓ Draining waste water during cycle.
- ✓ All purging system to remove air before sterile time.
- ✓ Total cycle time required only 12 mins.
- ✓ When inner pressure is over heated PLC will drain out Steam to reduce pressure.



Technical Specification:

Sterilizing temperature	121°C/134°C
Sterilization Time	4 mins
Sterilization Cycle Time	11 mins
Special sterilization	Kill active HV,HBV,BSE and spore
Dry procedure	Strong vacuum drying, the residual humidity<0.2%
Display	Widescreen LCD
Test	Bowie & Dick Steam penetration test, vacuum test , Helix test
Safety	Safety valves, Double Pressure lock Automatic check.
Water supply system	3L distiller water tank,3L waste water tank
Clean procedure	Automatic clean inside pipe and steam generator
Instrument containers	5 levels with 3 Round trays

STERIMAXX INC

"A reliable partner in Sterilization !!!

A/55 Shankesh Bldg, Sanghvi Nagar,
Mira Bhayander Road, Mira Road (East).
Thane 401 107 Maharashtra, INDIA.

Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com