

## Sterilization Flat Reel

### Material features:

Imported 60GSM Medical Grade paper certified BS 6256 & EN 868 Part – I & III & Tyvek I073B with heat stable special multilayered laminated film.

Can be easily sealed with all sealing machines.

Triple band seal creates three independent barriers to Contaminations, while reducing the risk of fiber tear.

Non-toxic accurate, process indicators for steam and EO gas conforming to the EN 868.

Strong Seals to eliminate bursting during autoclaving.



### Sterilization Flat Reels

Size In Inches	Size in Cm
2 inches	5 cm x 200 meters
3 inches	7.5 cm x 200 meters
4 inches	10 cm x 200 meters
6 inches	15 cm x 200 meters
8 inches	20 cm x 200 meters
10 inches	25 cm x 200 meters
12 inches	30 cm x 200 meters

## Self Seal Pouches

- Easy to use pouches, without a heat sealer.
- Pre folded flap for tight sealing.
- Steam & EO gas process indicators are imprints on pouches.

### Sterilization Self Seal Pouches

Size in mm
100 x 275mm
150 x 275mm
250 x 400 mm
330 x 450 mm



## Medical Pouch Sealing Machine

The machine is one of the matching products of the sterilizer with the features of stable performance, nice appearance, easy operation and low repair rate etc. The best sealing temperature can be adjusted as the voltage changes. Stainless steel is used for its shell to prevent corrosion and make it easy for cleaning. Cutter knife of double blades is installed inside and the cut is even. Seal strength of 10 mm.

### Technical Parameters

Source Voltage and Frequency: AC220V/50Hz,

Power: 500W

Greatest Sealing Width: 250mm

Indentation Width: 10mm Gross



**STERIMAXX INC**

*"A reliable partner in Sterilization !!!*

A/55 Shankesh Bldg, Sanghvi Nagar,  
Mira Bhayander Road, Mira Road (East).

Thane 401 107 Maharashtra, INDIA.

Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com