

## Ultrasonic Cleaner - ABS

- ✓ Easy to Operate.
- ✓ With Heater.
- ✓ 5 Recycle Digital timer for cleaning.
- ✓ Safety and reliable with Temp (65°C) & Time. (45mins) controller.
- ✓ Inbuilt Circuit Protector to increase unit life.
- ✓ Radiator fan for cooling to protect during long working hrs.
- ✓ Wet proof PCB.
- ✓ Strong Plastic Body.



### **Technical Specification:**

Rating voltage & Frequency	AC 220-240V ~ 50/60Hz
Power Consumption	170W
Ultrasonic Frequency	42,000 Hz
Material of Tank	Stainless Steel SUS304
Material of Plastic Housing	ABS 757
Tank Size (L x W x H)	264 x 164 x 80 mm
Tank Capacity	2600 ml
5 Recycle Digital Timer	90s/180s/280s/380s/480secs
Net Wt & Gross Wt	2.5Kg / 3.0 Kg



**STERIMAXX INC**

*"A reliable partner in Sterilization !!!*

A/55 Shankesh Bldg, Sanghvi Nagar,  
Mira Bhayander Road, Mira Road (East).  
Thane 401 107 Maharashtra, INDIA.  
Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com