

Ultrasonic Cleaner - ABS

- ✓ Control Panel with capacitive sensing technology.
- ✓ Durable, reliable, resistant to water and harsh chemicals.
- ✓ Industrial ultrasonic transducer (60W) with stronger cleaning effect.
- ✓ 2 color led display with 1 to 30 mins full range timer, special degas feature.
- ✓ Ceramic Heater of 100 Watt.
- ✓ Multiple circuit protectors to protect against improper handling and overload.
- ✓ Inbuilt cooling fan.
- ✓ Moisture proof PCB and industrial grade IC.
- ✓ Engineering grade plastic 4mm thick wall for better impact and water-proof performance.



Technical Specification:

Rating voltage & Frequency	AC 220-240V ~ 50/60Hz
Power Consumption	170W
Ultrasonic Frequency	35 KHz
Material of Tank	Stainless Steel SUS304
Material of Plastic Housing	ABS 757
Tank Size (L x W x H)	264 x 150 x 76 mm
Tank Capacity	2500 ml
Net Wt & Gross Wt	2.5Kg / 3.0 Kg



STERIMAXX INC

"A reliable partner in Sterilization !!!"

A/55 Shankesh Bldg, Sanghvi Nagar,
Mira Bhayander Road, Mira Road (East).
Thane 401 107 Maharashtra, INDIA.
Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com