



Ultrasonic Cleaner SS (LCD)



Multifunctional ultrasonic cleaner satisfies your cleaning requirement!

- With large LCD display to show the precise time, temperature, power level
- Time setting with 1-99 minutes and built-in heating to 60°C, adjustable power range of 10%~100%.
- Built-in "sound insulation paste" keep this series with low noise
- With sweep function & degas function
- Humanized carrying handle and drain valve for whole models
- Memory function will automatically save the most recently used setup parameters for next use.
- With sleeping mode to keep electrical energy saving



Adjustable power



Memory function



Low noise



Large LCD display



Sleeping mode



Sweep function



Adjustable time



Heating function



Degas function

Adjustable power

Power is adjustable (step by 10%) accordingly to the user's actual requirement



Large LCD display

Large LCD display screen always show the clear parameters when operate different functions



Sweep function

Make the cleaning process
is more evenly and more
thoroughly



Make the cleaning process
is more efficiently and
more quickly

This series is special



Built-in "sound insulation
paste" keep this series with low
noise

Specifications:

Model	Tank Capacity Liter	Tank Size in mm	Unit Size in mm	Ultrasonic Power(Max.) [power is adjustable]	Heating Power (Max.)
ATS-1	1.5	150x135x100	230x232x298	70W	60W
ATS-2	2.5	240x135x100	265x165x235	120W	100W
ATS-3	3.5	240x135x150	265x165x285	120W	100W
ATS-6	6.5	300x150x150	330x180x310	180W	200W
ATS-7	7	300x240x100	330x270x260	240W	500W
ATS-10	10	300x240x150	330x270x310	240W	500W
ATS-15	15	330x300x150	360x330x330	360W	500W
ATS-20	22	500x300x150	530x330x340	480W	500W
ATS-30	30	500x300x200	530x330x390	600W	500W
ATS-45	45	500x300x300	530x330x490	720W	1000W



Applications:

- Cleaning the glass ware
- Lab instruments
- Dissolve, disperse and emulsify samples
- Sonochemistry i.e. degassing the solution
- Degassing HPLC Samples
- Sonication

Why Our Ultrasonic Sonicator Bath?

- Use of Pulse Sweep Energy
- Precision cleaning
- Time saving
- Provide safe solution
- Highly effective for removing chemical residues from Ultrasonic Cleaner / Sonicator Bath



STERIMAXX INC

"A reliable partner in Sterilization !!!"

STERIMAXX INC

"A reliable partner in Sterilization !!!"

A/55 Shankesh Bldg, Sanghvi Nagar,
Mira Bhayander Road, Mira Road (East).
Thane 401 107 Maharashtra, INDIA.
Cell: +91 9987 51 4946.

sterimaxx@gmail.com. www.sterimax.com